

3Shape

Import guide



ZERAMEX

How to import the library into your 3Shape system?

1. Open the 3Shape Dental System Control Panel
2. Tools ► Import/Export ► Import material
3. Select “zeramex_xt_3shape_v5.dme”. Save and close the “3Shape Dental System Control Panel”.
4. You have successfully installed the ZerameX XT Zerabase library.

Please **contact the support service of your CAD/CAM (support@3shape.com) provider with any questions** on software data import.

Order the parts necessary for treatment directly with ZerameX.

Order-Hotline

T 00800 – 93 55 66 37

F 00800 – 93 55 63 77

order@zeramex.com

