

Zeramex XT Range

XT



ZERAMEX
naturally, white implants


| Zeramex XT Implants | | | | |
|---|----------|--|---------------|---------------------------|
| SB Ø3.5 mm | Art. No. | Name | Dimension | Material |
|  | XT15508 | Zeramex XT Implant Ø3.5×8 mm SB (incl. Healing Cap) | Length: 8 mm | ZrO ₂ -ATZ-HIP |
|  | XT15510 | Zeramex XT Implant Ø3.5×10 mm SB (incl. Healing Cap) | Length: 10 mm | |
|  | XT15512 | Zeramex XT Implant Ø3.5×12 mm SB (incl. Healing Cap) | Length: 12 mm | |
| RB Ø4.2 mm | Art. No. | Name | Dimension | Material |
|  | XT16508 | Zeramex XT Implant Ø4.2×8 mm RB (incl. Healing Cap) | Length: 8 mm | ZrO ₂ -ATZ-HIP |
|  | XT16510 | Zeramex XT Implant Ø4.2×10 mm RB (incl. Healing Cap) | Length: 10 mm | |
|  | XT16512 | Zeramex XT Implant Ø4.2×12 mm RB (incl. Healing Cap) | Length: 12 mm | |
|  | XT16514 | Zeramex XT Implant Ø4.2×14 mm RB (incl. Healing Cap) | Length: 14 mm | |
| WB Ø5.5 mm | Art. No. | Name | Dimension | Material |
|  | XT17508 | Zeramex XT Implant Ø5.5×8 mm WB (incl. Healing Cap) | Length: 8 mm | ZrO ₂ -ATZ-HIP |
|  | XT17510 | Zeramex XT Implant Ø5.5×10 mm WB (incl. Healing Cap) | Length: 10 mm | |
|  | XT17512 | Zeramex XT Implant Ø5.5×12 mm WB (incl. Healing Cap) | Length: 12 mm | |

Zeramex XT Abutments

| SB Platform | Art. No. | Name | Dimension | Material |
|-------------|----------|--|-----------------------------------|--------------------------------------|
| | SB15501 | Zeramex XT Abutment SB Straight, 1 mm (incl. screw) | AH: 7 mm CH: 1 mm Ø: 4.6 mm | ZrO ₂ -ATZ-HIP Vicarbo |
| | SB15502 | Zeramex XT Abutment SB Straight, 2 mm (incl. screw) | AH: 8 mm CH: 2 mm Ø: 4.6 mm | |
| | SB15515 | Zeramex XT Abutment SB Angular 15°, 1 mm (incl. screw) | AH: 7 mm CH: 1 mm Ø: 4.6 mm | |

| SB CADCAM | Art. No. | Name | Dimension | Material |
|-----------|----------|--|---------------------------------------|--------------------------------------|
| | SB15535 | Zeramex XT Zerabase X SB, Engaging for crowns (incl. screw) | AH: 4.8 mm CH: 0.8 mm Ø: 4.6 mm | ZrO ₂ -ATZ-HIP Vicarbo |
| | SB15536 | Zeramex XT Zerabase X SB Non-engaging for bars & bridges (incl. screw) | AH: 4.8 mm CH: 0.8 mm Ø: 4.6 mm | |

Zerabase Abutments are integrated into the systems of 3Shape and exocad.




| SB/RB/WB Vicarbo Screw | Art. No. | Name | Dimension | Material |
|------------------------|----------|--|----------------|----------|
| | RB16550 | Zeramex XT (SB/RB/WB) Vicarbo screw, Straight/Angular/CADCAM/Zerabase X/Customized abutments | Length: 7.4 mm | Vicarbo |

Zeramex XT Abutments

| RB Platform | Art. No. | Name | Dimension | Material |
|-------------|----------|--|---------------------------------|--------------------------------------|
| | RB16501 | Zeramex XT Abutment RB Straight, 1 mm (incl. screw) | AH: 7 mm CH: 1 mm Ø: 5 mm | ZrO ₂ -ATZ-HIP Vicarbo |
| | RB16502 | Zeramex XT Abutment RB Straight, 2 mm (incl. screw) | AH: 8 mm CH: 2 mm Ø: 5 mm | |
| | RB16515 | Zeramex XT Abutment RB Angular 15°, 1 mm (incl. screw) | AH: 7 mm CH: 1 mm Ø: 5 mm | |

| RB CADCAM | Art. No. | Name | Dimension | Material |
|-----------|----------|--|---------------------------------------|--------------------------------------|
| | RB16535 | Zeramex XT Zerabase X RB, Engaging for crowns (incl. screw) | AH: 4.8 mm CH: 0.8 mm Ø: 4.6 mm | ZrO ₂ -ATZ-HIP Vicarbo |
| | RB16536 | Zeramex XT Zerabase X RB Non-engaging for bars & bridges (incl. screw) | AH: 4.8 mm CH: 0.8 mm Ø: 4.6 mm | |
| | RB16530 | Zeramex XT Zerabase RB, Engaging for crowns (incl. screw) | AH: 7 mm CH: 1 mm Ø: 5 mm | |
| | RB16531 | Zeramex XT Zerabase RB Non-engaging for bars & bridges (incl. screw) | AH: 7 mm CH: 1 mm Ø: 5 mm | |

Zerabase X Abutments are integrated into the systems of 3Shape and exocad.










| SB/RB/WB Vicarbo Screw | Art. No. | Name | Dimension | Material |
|------------------------|----------|--|----------------|----------|
| | RB16550 | Zeramex XT (SB/RB/WB) Vicarbo screw, Straight/Angular/CADCAM/Zerabase X/Customized abutments | Length: 7.4 mm | Vicarbo |

Range


XT


Zeramex XT Abutments

| WB Platform | Art. No. | Name | Dimension | Material |
|---|----------|--|---------------------------------|--------------------------------------|
|  | WB17501 | Zeramex XT Abutment WB Straight, 1 mm (incl. screw) | AH: 7 mm CH: 1 mm Ø: 6 mm | ZrO ₂ -ATZ-HIP Vicarbo |
|  | WB17502 | Zeramex XT Abutment WB Straight, 2 mm (incl. screw) | AH: 8 mm CH: 2 mm Ø: 6 mm | |
|  | WB17515 | Zeramex XT Abutment WB Angular 15°, 1 mm (incl. screw) | AH: 7 mm CH: 1 mm Ø: 6 mm | |




| WB CADCAM | Art. No. | Name | Dimension | Material |
|---|----------|--|---------------------------------------|--------------------------------------|
|  | WB17535 | Zeramex XT Zerabase X WB Engaging for crowns (incl. screw) | AH: 4.8 mm CH: 0.8 mm Ø: 5.6 mm | ZrO ₂ -ATZ-HIP Vicarbo |
|  | WB17536 | Zeramex XT Zerabase X WB Non-engaging for bars & bridges (incl. screw) | AH: 4.8 mm CH: 0.8 mm Ø: 5.6 mm | |
|  | WB17530 | Zeramex XT Zerabase WB Engaging for crowns (incl. screw) | AH: 7 mm CH: 1 mm Ø: 6 mm | |
|  | WB17531 | Zeramex XT Zerabase WB Non-engaging for bars & bridges (incl. screw) | AH: 7 mm CH: 1 mm Ø: 6 mm | |






Zerabase Abutments are integrated into the systems of 3Shape and exocad.



| SB/RB/WB Vicarbo Screw | Art. No. | Name | Dimension | Material |
|---|----------|--|----------------|----------|
|  | RB16550 | Zeramex XT (SB/RB/WB) Vicarbo screw, Straight/Angular/CADCAM/Zerabase X/Customized abutments | Length: 7.4 mm | Vicarbo |

Zeramex Prosthetics


| SB soft tissue management | Art. No. | Name | Dimension | Material |
|---|----------|--|---------------------------------|-----------------|
|  | SB35500 | Zeramex XT Healing Cap SB | Height: 0.6 mm | PEEK |
|  | SB35503 | Zeramex XT Gingivaformer SB, 3 mm | Height: 3 mm Ø: 4 mm | |
|  | SB35504 | Zeramex XT Gingivaformer SB, 4 mm | Height: 4 mm Ø: 4 mm | |
|  | SB35530 | Zeramex XT Provisional SB (incl. Provisional Screw) | AH: 7 mm CH: 1 mm Ø: 4 mm | PEEK Vicarbo |






| SB Impression Taking | Art. No. | Name | Dimension | Material |
|---|----------|---|--|------------------------|
|  | SB35510 | Zeramex XT Transfer Open Tray SB | High sleeve: 11 mm / height incl. screw: 20 mm | Aluminium PEEK-CW30 |
|  | SB35512 | Zeramex XT Transfer Closed Tray SB | High sleeve: 7 mm / height incl. screw: 14 mm | |
|  | SB35513 | Zeramex XT Transfer Closed Tray, long SB | High sleeve: 11 mm / height incl. screw: 18 mm | |
|  | SB35514 | Zeramex XT Scanbody SB (incl. Provisional Screw) | Height: 10 mm | PEEK Vicarbo |
|  | SB35522 | Zeramex XT Digital Implant Replica SB | Length: 10 mm | Aluminium |




| RB soft tissue management | Art. No. | Name | Dimension | Material |
|---|----------|-----------------------------------|-------------------------|----------|
|  | RB36500 | Zeramex XT Healing Cap RB | Height: 0.6 mm | PEEK |
|  | RB36503 | Zeramex XT Gingivaformer RB, 3 mm | Height: 3 mm Ø: 5 mm | |
|  | RB36504 | Zeramex XT Gingivaformer RB, 4 mm | Height: 4 mm Ø: 5 mm | |




Range

XT

| | | | | |
|---|---------|--|---------------------------------|-----------------|
|  | RB36530 | Zeramex XT Provisional RB (incl. Provisional Screw) | AH: 7 mm CH: 1 mm Ø: 5 mm | PEEK Vicarbo |
|---|---------|--|---------------------------------|-----------------|



| RB Impression Taking | Art. No. | Name | Dimension | Material |
|--|----------|---|--|------------------------|
|  | RB36510 | Zeramex XT Transfer Open Tray RB | High sleeve: 11 mm / height incl. screw: 20 mm | Aluminium PEEK-CW30 |
|  | RB36512 | Zeramex XT Transfer Closed Tray RB | High sleeve: 7 mm / height incl. screw: 14 mm | |
|  | RB36513 | Zeramex XT Transfer Closed Tray, long RB | High sleeve: 11 mm / height incl. screw: 18 mm | |
|  | RB36514 | Zeramex XT Scanbody RB (incl. Provisional Screw) | Height: 10 mm | PEEK Vicarbo |
|  | RB36522 | Zeramex XT Digital Implant Replica RB | Length: 10 mm | Aluminium |




| WB soft tissue management | Art. No. | Name | Dimension | Material |
|---|----------|--|---------------------------------|-----------------|
|  | WB37500 | Zeramex XT Healing Cap WB | Height: 0.6 mm | PEEK |
|  | WB37503 | Zeramex XT Gingivaformer WB, 3 mm | Height: 3 mm Ø: 6 mm | |
|  | WB37504 | Zeramex XT Gingivaformer WB, 4 mm | Height: 4 mm Ø: 6 mm | |
|  | WB37530 | Zeramex XT Provisional WB (incl. Provisional Screw) | AH: 7 mm CH: 1 mm Ø: 6 mm | PEEK Vicarbo |





| WB Impression Taking | Art. No. | Name | Dimension | Material |
|---|----------|--|--|------------------------|
|  | WB37510 | Zeramex XT Transfer Open Tray WB | High sleeve: 11 mm / height incl. screw: 20 mm | Aluminium PEEK-CW30 |
|  | WB37512 | Zeramex XT Transfer Closed Tray WB | High sleeve: 7 mm / height incl. screw: 14 mm | |
|  | WB37513 | Zeramex XT Transfer Closed Tray, long WB | High sleeve: 11 mm / height incl. screw: 18 mm | |

Range

XT







| | | | | |
|---|---------|---|---------------|-------------------|
|  | WB37514 | Zeramex XT Scanbody WB (incl. screw) | Height: 10 mm | PEEK PEEK-CW30 |
|  | WB37522 | Zeramex XT Digital Implant Replica WB | Length: 10 mm | Aluminium |

| SB/RB/WB Auxiliary Parts | Art. No. | Name | Dimension | Material |
|---|----------|---|----------------|-----------|
|  | RB36521 | Zeramex XT Digital Implant Replica Placer (SB/RB/WB) | Length: 14 mm | PEEK-CW30 |
|  | RB36550 | Zeramex XT Provisional Screw (SB/RB/WB) | Length: 8.6 mm | Vicarbo |
|  | RB36554 | Zeramex XT Try-In Screw (SB/RB/WB) | Length: 7.1 mm | PEEK |

| Zeramex Docklocs® | Art. No. | Name | Material |
|---|----------|---|--|
|  | SB15542 | Zeramex Docklocs® Abutment, 2 mm | ZrO ₂ -ATZ-HIP Vicarbo |
|  | SB15543 | Zeramex Docklocs® Abutment, 3 mm | |
|  | SB15544 | Zeramex Docklocs® Abutment, 4 mm | |
|  | XT38227 | Zeramex Docklocs® Insertion Instrument | Stainless steel PEEK |
|  | XT38253 | Docklocs® Laboratory Set, up to 40° divergence compensation: 2 pcs zirconia oxide retention housing (Ø5.8 mm, height 2.5 mm) with black processing insert (height 1.9 mm), 2 pcs block-out ring, 2 pcs replacement male, green, 2 pcs replacement male, orange, 2 pcs replacement male, red | Santroprene® TPE Polyamide Housing ZrO ₂ HD-PE Purell |
|  | XT38251 | Docklocs® Laboratory Set, up to 40° divergence compensation: 2 pcs titanium retention housing (Ø5.5 mm, height 2.5 mm) with black processing insert (height 1.9 mm), 2 pcs block-out ring, 2 pcs replacement male, green, 2 pcs replacement male, orange, 2 pcs replacement male, red | Santroprene® TPE Polyamide Titanium housing HD-PE Purell |
















Range

| | | | |
|---|---------|--|--|
|  | XT38205 | Docklocs® Replacement Male, red Extra-light retention, 10°–20° 8 pcs | Polyamide |
|  | XT38206 | Docklocs® Replacement Male, orange Light retention, 10°–20° 8 pcs | |
|  | XT38207 | Docklocs® Replacement Male, green Strong retention, 10°–20° 8 pcs | |
|  | XT38209 | Docklocs® Block-out Ring 20 pcs | Santoprene® TPE |
|  | XT38230 | Zirconia oxide retention housing with processing insert 2 pcs | Housing ZrO ₂ HD-PE Purell |
|  | XT38210 | Titanium retention housing with processing insert 4 pcs | Titanium housing G5 HD-PE Purell |
|  | XT38214 | Docklocs® lab analogue straight (Ø4 mm) 4 pcs | Grade 5 titanium |
|  | XT38215 | Docklocs® impression coping with black processing insert 4 pcs | Titanium housing G5 HD-PE Purell |

| SB/RB/WB Tools | Art. No. | Name | Material |
|---|----------|-----------------------------|--|
|  | P35601 | Zeramex Rosedrill Ø2 mm | Stainless steel |
|  | T35602 | Zeradrill Pilot Ø2.3 mm | |
|  | XT35630 | Zeramex Profile Drill SB | Stainless steel with carbon coating |
|  | T35608 | Zeradrill S8 (Small 8 mm) | |
|  | T35610 | Zeradrill S10 (Small 10 mm) | |
|  | T35612 | Zeradrill S12 (Small 12 mm) | |















Range



| | | | |
|---|--------|-------------------------------|-------------------------------------|
|  | T35614 | Zeradrill S14 (Small 14 mm) | Stainless steel with carbon coating |
|  | T35620 | Zeratap Small Ø3.5 mm | |
|  | T35622 | Zeradrill Extension Small | |
|  | T36608 | Zeradrill R8 (Regular 8 mm) | |
|  | T36610 | Zeradrill R10 (Regular 10 mm) | |
|  | T36612 | Zeradrill R12 (Regular 12 mm) | |
|  | T36614 | Zeradrill R14 (Regular 14 mm) | |
|  | T36620 | Zeratap Regular Ø4.2 mm | |
|  | T36622 | Zeradrill Extension Regular | |
|  | T37608 | Zeradrill W8 (Wide 8 mm) | |
|  | T37610 | Zeradrill W10 (Wide 10 mm) | |
|  | T37612 | Zeradrill W12 (Wide 12 mm) | |
|  | T37620 | Zeratap Wide Ø5.5 mm | |
|  | T37622 | Zeradrill Extension Wide | |
|  | C7650* | Drill extension | |


*Art. KI589B for Switzerland












Range

| | | | |
|---|---------|--|-----------------|
|  | T38650 | Zeramex T depth gauge, 4 pcs | Stainless steel |
|  | XT38619 | Zeramex Prosthetic Key, 19 mm | |
|  | XT38623 | Zeramex Prosthetic Key, 23 mm | |
|  | XT38628 | Zeramex Prosthetic Key, 28 mm | |
|  | XT36620 | Zeramex Pick-up, 20 mm SB/RB/WB | |
|  | XT36625 | Zeramex Pick-up, 25 mm SB/RB/WB | |
|  | XT36622 | Zeramex Rescue Pick-up SB/RB/WB | |
|  | P48932 | Zeramex P Ratchet Adapter Unit Short | |
|  | P48935 | Zeramex P Surgical Ratchet without Adapter | |
|  | XT35651 | Zeramex XT Drill Stop, for Small Drill | |
|  | XT36651 | Zeramex XT Drill Stop, for Regular Drill | |
|  | XT37651 | Zeramex XT Drill Stop, for Wide Drill | |
|  | XT48860 | Zeramex XT Prosthetics Kit | - |
|  | XT48865 | Zeramex XT Prosthetics Kit incl. Ratchet | |















Range



| | | | |
|---|---------|--|---|
|  | XT48850 | Zeramex XT Surgery Tray | |
| | XT48854 | Zeramex XT Surgery Tray, fully equipped, incl. Ratchet | - |

| Zeramex Digital Solutions* | Art. No. | Name | Material |
|---|----------|--|----------------------------------|
|  | SB15551 | Customized Abutment SB (incl. Vicarbo screw) | ZrO ₂ -TZP Vicarbo |
|  | SB15552 | Customized zirconia crown for Zerabase <u>without</u> screw hole (unglazed & excl. Zerabase Abutment) SB | ZrO ₂ -TZP |
|  | SB15553 | Customized zirconia crown for Zerabase <u>with</u> screw hole (unglazed & excl. Zerabase Abutment) SB | |
|  | SB15554 | Customized zirconia Gingivaformer SB (incl. Vicarbo screw) | ZrO ₂ -TZP Vicarbo |
|  | SB15560 | Monolithic crown SB (incl. Vicarbo screw) | ZrO ₂ -TZP Vicarbo |
|  | SB15561 | 2 monolithic crowns splinted SB (incl. Vicarbo screw) | |
|  | SB15565 | Monolithic bridge 3-unit SB (incl. Vicarbo screw) | |
|  | SB15566 | Monolithic bridge 3-unit SB (incl. Vicarbo screw) cutback | |
|  | SB15590 | Package customized: Zirconia Gingivaformer & Abutment SB (incl. Vicarbo screw) | |
|  | SB15591 | Package Monolithic: Zirconia Gingivaformer & crown SB (incl. Vicarbo screw) | |
|  | RB16551 | Customized Abutment RB (incl. Vicarbo screw) | |










Range

| | | | |
|---|---------|--|----------------------------------|
|  | RB16552 | Customized zirconia crown for Zerabase <u>without</u> screw hole (unglazed & excl. Zerabase Abutment) RB | ZrO ₂ -TZP |
|  | RB16553 | Customized zirconia crown for Zerabase <u>with</u> screw hole (unglazed & excl. Zerabase Abutment) RB | |
|  | RB16554 | Customized zirconia Gingivaformer RB (incl. Vicarbo screw) | ZrO ₂ -TZP Vicarbo |
|  | RB16560 | Monolithic crown RB (incl. Vicarbo screw) | ZrO ₂ -TZP Vicarbo |
|  | RB16561 | 2 monolithic crowns splinted RB (incl. Vicarbo screw) | |
|  | RB16565 | Monolithic bridge 3-unit RB (incl. Vicarbo screw) | |
|  | RB16566 | Monolithic bridge 3-unit RB (incl. Vicarbo screw) cutback | |
|  | RB16590 | Package customized: Zirconia Gingivaformer & Abutment RB (incl. Vicarbo screw) | |
|  | RB16591 | Package Monolithic: Zirconia Gingivaformer & crown RB (incl. Vicarbo screw) | |
|  | WB17551 | Customized Abutment WB (incl. Vicarbo screw) | |
|  | WB17552 | Customized zirconia crown for Zerabase <u>without</u> screw hole (unglazed & excl. Zerabase Abutment) WB | ZrO ₂ -TZP |
|  | WB17553 | Customized zirconia crown for Zerabase <u>with</u> screw hole (unglazed & excl. Zerabase Abutment) WB | |
|  | WB17554 | Customized zirconia Gingivaformer WB (incl. Vicarbo screw) | ZrO ₂ -TZP Vicarbo |
|  | WB17560 | Monolithic crown WB (incl. Vicarbo screw) | ZrO ₂ -TZP Vicarbo |
|  | WB17561 | 2 monolithic crowns splinted WB (incl. Vicarbo screw) | |

Range



| | | | |
|---|---------|--|----------------------------------|
|  | WB17565 | Monolithic bridge 3-unit WB (incl. Vicarbo screw) | ZrO ₂ -TZP Vicarbo |
|  | WB17566 | Monolithic bridge 3-unit WB (incl. Vicarbo screw) cutback | |
|  | WB17590 | Package customized: Zirconia Gingivaformer & Abutment WB (incl. Vicarbo screw) | |
|  | WB17591 | Package Monolithic: Zirconia Gingivaformer & crown WB (incl. Vicarbo screw) | |







| Zeramex Digital Solutions Multilayer* | Art. No. | Name | Material |
|---|----------|--|----------------------------------|
|  | SB15751 | Customized Abutment SB, ML (incl. Vicarbo screw) | ZrO ₂ -TZP Vicarbo |
|  | SB15752 | Customized zirconia crown for Zerabase <u>without</u> screw hole, ML (unglazed & excl. Zerabase Abutment) SB | ZrO ₂ -TZP |
|  | SB15753 | Customized zirconia crown for Zerabase <u>with</u> screw hole, ML (unglazed & excl. Zerabase Abutment) SB | |
|  | SB15754 | Customized zirconia Gingivaformer SB, ML (incl. Vicarbo screw) | ZrO ₂ -TZP Vicarbo |
|  | SB15760 | Monolithic crown SB, ML (incl. Vicarbo screw) | ZrO ₂ -TZP Vicarbo |
|  | SB15761 | 2 monolithic crowns splinted SB, ML (incl. Vicarbo screw) | |
|  | SB15765 | Monolithic bridge 3-unit SB, ML (incl. Vicarbo screw) | |
|  | SB15766 | Monolithic bridge 3-unit SB, ML (incl. Vicarbo screw) cutback | |
|  | SB15790 | Package customized: Zirconia Gingivaformer & Abutment SB, ML (incl. Vicarbo screw) | |



Range

| | | | |
|---|---------|--|----------------------------------|
|  | SB15791 | Package Monolithic: Zirconia Gingivaformer & crown SB, ML (incl. Vicarbo screw) | ZrO ₂ -TZP Vicarbo |
|  | RB16751 | Customized Abutment RB, ML (incl. Vicarbo screw) | |
|  | RB16752 | Customized zirconia crown for Zerabase <u>without</u> screw hole, ML (unglazed & excl. Zerabase Abutment) RB | ZrO ₂ -TZP |
|  | RB16753 | Customized zirconia crown for Zerabase <u>with</u> screw hole, ML (unglazed & excl. Zerabase Abutment) RB | |
|  | RB16754 | Customized zirconia Gingivaformer RB, ML (incl. Vicarbo screw) | ZrO ₂ -TZP Vicarbo |
|  | RB16760 | Monolithic crown RB, ML (incl. Vicarbo screw) | ZrO ₂ -TZP Vicarbo |
|  | RB16761 | 2 monolithic crowns splinted RB, ML (incl. Vicarbo screw) | |
|  | RB16765 | Monolithic bridge 3-unit RB, ML (incl. Vicarbo screw) | |
|  | RB16766 | Monolithic bridge 3-unit RB, ML (incl. Vicarbo screw) cutback | |
|  | RB16790 | Package customized: Zirconia Gingivaformer & Abutment RB, ML (incl. Vicarbo screw) | |
|  | RB16791 | Package Monolithic: Zirconia Gingivaformer & crown RB, ML (incl. Vicarbo screw) | |
|  | WB17751 | Customized Abutment WB, ML (incl. Vicarbo screw) | ZrO ₂ -TZP |
|  | WB17752 | Customized zirconia crown for Zerabase <u>without</u> screw hole, ML (unglazed & excl. Zerabase Abutment) WB | |
|  | WB17753 | Customized zirconia crown for Zerabase <u>with</u> screw hole, ML (unglazed & excl. Zerabase Abutment) WB | |
|  | WB17754 | Customized zirconia Gingivaformer WB, ML (incl. Vicarbo screw) | ZrO ₂ -TZP Vicarbo |

Range



| | | | |
|---|---------|---|----------------------------------|
|  | WB17760 | Monolithic crown WB, ML (incl. Vicarbo screw) | ZrO ₂ -TZP Vicarbo |
|  | WB17761 | 2 monolithic crowns splinted WB, ML (incl. Vicarbo screw) | |
|  | WB17765 | Monolithic bridge 3-unit WB, ML (incl. Vicarbo screw) | |
|  | WB17766 | Monolithic bridge 3-unit WB, ML (incl. Vicarbo screw) cutback | |
|  | WB17790 | Package customized: Zirconia Gingivaformer & Abutment WB, ML (incl. Vicarbo screw) | |
|  | WB17791 | Package Monolithic: Zirconia Gingivaformer & crown WB, ML (incl. Vicarbo screw) | |

| Zeramex Digital Solutions Services* | Art. No. | Name | Material |
|---|----------|--|----------|
|  | RB16570 | 3D model OK/UK (without Replica) SB/RB/WB | |
|  | RB16571 | Stain & Glaze SB/RB/WB | |
| - | RB16580 | Zeramex SB, RB & WB Individuelle Designentwicklung | - |
| - | RB16581 | Zeramex Scan (alle Produkte) | |
| - | RB16582 | Zeramex SB, RB & WB Virtuelles Modelldesign aufgrund Intraoralscandaten | |

***Sample images/differs from original!**

Docklocs® is a registered trademark of MEDEALIS GmbH, DE.
Santoprene® is a registered trademark of Exxon Mobil Corporation, USA.

General Information

Guarantee

CeramTec Switzerland LLC offers a lifelong guarantee for implants, and a 10-year guarantee for abutments and Vicarbo screws. The Zeramex warranty conditions apply.

Certification

CE / CE₀₄₈₃ / CE₀₅₀₀ (see p. 68)

Delivery and packaging

Delivery is in accordance with the general terms and conditions (T&Cs) of CeramTec Switzerland LLC. Intact sterile packaging protects the implant from external influences and ensures sterile storage up to the printed expiration date. Zeramex XT implants and components must be stored dry in their original packaging at room temperature and protected from sunlight. Only open the packaging shortly before surgery. We recommend comprehensive clinical, radiological and statistical documentation. The inside labels (patient label) must allow traceability of the implants.

Exclusion of liability

Zeramex XT implants are part of an overall system and may be used only with the components designed

for this system. CeramTec Switzerland LLC will not be held liable for any damage arising from improper use, or from using non-original components. The general terms and conditions of CeramTec Switzerland LLC also apply.

Training

For information on courses and further education for the Zeramex XT System, please contact us at www.zeramex.com.

Material properties

All implants and abutments are made from hot-densified zirconia oxide ATZ-HIP® (HIP = Hot Isostatic Postcompaction). For reasons of quality and strength, the implants and abutments are strictly machined into their final shape from solid, hard blanks using diamond-coated tools. The workpiece does not need finishing. This allows for highly precise and reproducible production of implants and abutments with the necessary precise fit.

Zerafil implant surface

- Microstructured
- Blasted and etched
- Hydrophilic

ZrO₂ ATZ-HIP

Zirconia dioxide, ATZ (alumina-toughened zirconia) (radiopaque)

Composition:

ZrO₂ 76%, Al₂O₃ 20%, Y₂O₃ 4%

Flexural strength: >1700 MPa (on average 2000 MPa)

CTE: 9x10⁻⁶/K

ZrO₂ TZP-A

Zirconia dioxide, TZP (tetragonal zirconia polycrystal) (radiopaque)

Composition:

ZrO₂ 95%, Al₂O₃ 5%, Y₂O₃ 0.25%

Flexural strength: 1,200 MPa.

CTE: 10,5

PEEK CLASSIX

Polyether ether ketone USP Class VI (not radiopaque)

Aluminium

Aluminium (not radiopaque)

PEEK CLASSIX CW30 LSG

Short carbon fibres (CF) in a PEEK CLASSIX LSG matrix (not radiopaque)

Composition:

CF 30%, PEEK Classix LSG 70%

Flexural strength: >130 MPa

Vicarbo

Unidirectional carbon fibres (CF) in a PEEK matrix (not radiopaque)

Composition:

CF 60%, PEEK 40%

Flexural strength: >900 MPa.

General Information



| | | | |
|---------|---|--|---|
| | CE Mark for MDR (2017/745) Class I devices | | CE Mark for legacy devices Class IIa / IIb |
| | CE Mark for MDR (2017/745) Class IIa / IIb devices | | Caution |
| | Article number | | Expiry date |
| | Batch code | | Observe the instructions for use |
| | Do not reuse | | Date of manufacture |
| | Non-sterile | | Manufacturer |
| | Sterilisation by steam or dry heat | | Do not use if the packaging is damaged |
| | Keep away from direct sunlight | | |
| | Keep dry | | Medical device |
| | Unique Device Identifier | | EU authorised representative |
| | Sterile barrier system | | Double packaging with single internal sterile barrier |
| | Do not resterilise | | Gamma sterilisation |
| Rx only | CAUTION: US federal law restricts the sale of this product to a physician or on his or her order. | | |
| | Implant Small ϕ 3.5 mm SB | | Implant Regular ϕ 4.2 mm RB |
| | Implant Wide ϕ 5.5 mm WB | | |

Contact

Europe

Order hotline for Germany, Switzerland and Austria

Phone: 00800 93 55 66 37

Fax: 00800 93 55 63 77

Email: order@zeramex.com

Zeramex Digital Solutions

Phone: 00800 04 00 13 33

Email: digitalsolutions@zeramex.com

Switzerland

CeramTec Switzerland LLC

Bodenackerstrasse 5

8957 Spreitenbach/Switzerland

Phone: 0041 44 388 36 36

Fax: 0041 44 388 36 39

Email: info@zeramex.com

Web: www.zeramex.com

Germany (EU sales)

CeramTec Dentalvertriebs GmbH

Wallbrunnstrasse 24

79539 Lörrach, Germany

Phone: +49 7621 1612 749

Fax: +49 7621 1612 780

E-mail: info@zeramex.com

Web: www.zeramex.com

Global

You can find the contact details of all sales partners on our website www.zeramex.com/contact



Zeramex XT

For natural red-white aesthetics





