

Zeramex XT Range

XT



ZERAMEX
naturally, white implants

Range

XT




Zeramex XT Implants


SB Ø3.5 mm	Art. No.	Name	Dimension	Material
	XT15508	Zeramex XT Implant Ø3.5×8 mm SB (incl. Healing Cap)	Length: 8 mm	ZrO ₂ -ATZ-HIP
	XT15510	Zeramex XT Implant Ø3.5×10 mm SB (incl. Healing Cap)	Length: 10 mm	
	XT15512	Zeramex XT Implant Ø3.5×12 mm SB (incl. Healing Cap)	Length: 12 mm	
RB Ø4.2 mm	Art. No.	Name	Dimension	Material
	XT16508	Zeramex XT Implant Ø4.2×8 mm RB (incl. Healing Cap)	Length: 8 mm	ZrO ₂ -ATZ-HIP
	XT16510	Zeramex XT Implant Ø4.2×10 mm RB (incl. Healing Cap)	Length: 10 mm	
	XT16512	Zeramex XT Implant Ø4.2×12 mm RB (incl. Healing Cap)	Length: 12 mm	
	XT16514	Zeramex XT Implant Ø4.2×14 mm RB (incl. Healing Cap)	Length: 14 mm	
WB Ø5.5 mm	Art. No.	Name	Dimension	Material
	XT17508	Zeramex XT Implant Ø5.5×8 mm WB (incl. Healing Cap)	Length: 8 mm	ZrO ₂ -ATZ-HIP
	XT17510	Zeramex XT Implant Ø5.5×10 mm WB (incl. Healing Cap)	Length: 10 mm	
	XT17512	Zeramex XT Implant Ø5.5×12 mm WB (incl. Healing Cap)	Length: 12 mm	

Range


XT

Zeramex XT Abutments

SB Platform	Art. No.	Name	Dimension	Material
	SB15501	Zeramex XT Abutment SB Straight, 1 mm (incl. screw)	AH: 7 mm CH: 1 mm Ø: 4.6 mm	ZrO ₂ -ATZ-HIP Vicarbo
	SB15502	Zeramex XT Abutment SB Straight, 2 mm (incl. screw)	AH: 8 mm CH: 2 mm Ø: 4.6 mm	
	SB15515	Zeramex XT Abutment SB Angular 15°, 1 mm (incl. screw)	AH: 7 mm CH: 1 mm Ø: 4.6 mm	






AH: Abutment height
CH: Collar height
Ø: Diameter


SB/RB/WB Vicarbo Screw	Art. No.	Name	Dimension	Material
	RB16550	Zeramex XT (SB/RB/WB) Vicarbo screw, Straight/Angular	Length: 7.4 mm	Vicarbo

Range


XT

Zeramex XT Abutments

RB Platform	Art. No.	Name	Dimension	Material
	RB16501	Zeramex XT Abutment RB Straight, 1 mm (incl. screw)	AH: 7 mm CH: 1 mm Ø: 5 mm	ZrO ₂ -ATZ-HIP Vicarbo
	RB16502	Zeramex XT Abutment RB Straight, 2 mm (incl. screw)	AH: 8 mm CH: 2 mm Ø: 5 mm	
	RB16515	Zeramex XT Abutment RB Angular 15°, 1 mm (incl. screw)	AH: 7 mm CH: 1 mm Ø: 5 mm	



AH: Abutment height
CH: Collar height
Ø: Diameter


SB/RB/WB Vicarbo Screw	Art. No.	Name	Dimension	Material
	RB16550	Zeramex XT (SB/RB/WB) Vicarbo screw, Straight/Angular	Length: 7.4 mm	Vicarbo

Range

XT

Zeramex XT Abutments

WB Platform	Art. No.	Name	Dimension	Material
	WB17501	Zeramex XT Abutment WB Straight, 1 mm (incl. screw)	AH: 7 mm CH: 1 mm Ø: 6 mm	ZrO ₂ -ATZ-HIP Vicarbo
	WB17502	Zeramex XT Abutment WB Straight, 2 mm (incl. screw)	AH: 8 mm CH: 2 mm Ø: 6 mm	
	WB17515	Zeramex XT Abutment WB Angular 15°, 1 mm (incl. screw)	AH: 7 mm CH: 1 mm Ø: 6 mm	



AH: Abutment height
CH: Collar height
Ø: Diameter

SB/RB/WB Vicarbo Screw	Art. No.	Name	Dimension	Material
	RB16550	Zeramex XT (SB/RB/WB) Vicarbo screw, Straight/Angular	Length: 7.4 mm	Vicarbo

Range

XT

Zeramex Prosthetics

SB soft tissue management	Art. No.	Name	Dimension	Material
	SB35500	Zeramex XT Healing Cap SB	Height: 0.6 mm	PEEK
	SB35503	Zeramex XT Gingivaformer SB, 3 mm	Height: 3 mm Ø: 4 mm	
	SB35504	Zeramex XT Gingivaformer SB, 4 mm	Height: 4 mm Ø: 4 mm	
	SB35530	Zeramex XT Provisional SB (incl. screw)	AH: 7 mm CH: 1 mm Ø: 4 mm	PEEK Vicarbo





SB Impression Taking	Art. No.	Name	Dimension	Material
	SB35510	Zeramex XT Transfer Open Tray SB	High sleeve: 11 mm / height incl. screw: 20 mm	Aluminium PEEK-CW30
	SB35512	Zeramex XT Transfer Closed Tray SB	High sleeve: 7 mm / height incl. screw: 14 mm	
	SB35513	Zeramex XT Transfer Closed Tray, long SB	High sleeve: 11 mm / height incl. screw: 18 mm	
	SB35522	Zeramex XT Implant Replica SB	Length: 10 mm	Aluminium

RB soft tissue management	Art. No.	Name	Dimension	Material
	RB36500	Zeramex XT Healing Cap RB	Height: 0.6 mm	PEEK
	RB36503	Zeramex XT Gingivaformer RB, 3 mm	Height: 3 mm Ø: 5 mm	
	RB36504	Zeramex XT Gingivaformer RB, 4 mm	Height: 4 mm Ø: 5 mm	




Range

XT

	RB36530	Zeramex XT Provisional RB (incl. screw)	AH: 7 mm CH: 1 mm Ø: 5 mm	PEEK Vicarbo
---	---------	---	---------------------------------	-----------------

RB Impression Taking	Art. No.	Name	Dimension	Material
	RB36510	Zeramex XT Transfer Open Tray RB	High sleeve: 11 mm / height incl. screw: 20 mm	Aluminium PEEK-CW30
	RB36512	Zeramex XT Transfer Closed Tray RB	High sleeve: 7 mm / height incl. screw: 14 mm	
	RB36513	Zeramex XT Transfer Closed Tray, long RB	High sleeve: 11 mm / height incl. screw: 18 mm	
	RB36522	Zeramex XT Implant Replica RB	Length: 10 mm	Aluminium




WB soft tissue management	Art. No.	Name	Dimension	Material
	WB37500	Zeramex XT Healing Cap WB	Height: 0.6 mm	PEEK
	WB37503	Zeramex XT Gingivaformer WB, 3 mm	Height: 3 mm Ø: 6 mm	
	WB37504	Zeramex XT Gingivaformer WB, 4 mm	Height: 4 mm Ø: 6 mm	
	WB37530	Zeramex XT Provisional WB (incl. screw)	AH: 7 mm CH: 1 mm Ø: 6 mm	PEEK Vicarbo







WB Impression Taking	Art. No.	Name	Dimension	Material
	WB37510	Zeramex XT Transfer Open Tray WB	High sleeve: 11 mm / height incl. screw: 20 mm	Aluminium PEEK-CW30
	WB37512	Zeramex XT Transfer Closed Tray WB	High sleeve: 7 mm / height incl. screw: 14 mm	
	WB37513	Zeramex XT Transfer Closed Tray, long WB	High sleeve: 11 mm / height incl. screw: 18 mm	

Range

XT

	WB37522	Zeramex XT Implant Replica WB	Length: 10 mm	Aluminium
---	---------	-------------------------------	---------------	-----------







SB/RB/WB Auxiliary Parts	Art. No.	Name	Dimension	Material
	RB36521	Zeramex XT Implant Replica Placer (SB/RB/WB)	Length: 14 mm	PEEK-CW30
	RB36550	Zeramex XT Provisional Screw (SB/RB/WB)	Length: 8.6 mm	Vicarbo
	RB36553	Zeramex XT Lab Screw (SB/RB/WB)	Length: 7.4 mm	PEEK

Zeramex Docklocs®	Art. No.	Name	Material
	SB15542	Zeramex Docklocs® Abutment, 2 mm	ZrO ₂ -ATZ-HIP Vicarbo
	SB15543	Zeramex Docklocs® Abutment, 3 mm	
	SB15544	Zeramex Docklocs® Abutment, 4 mm	
	XT38227	Zeramex Docklocs® Insertion Instrument	Stainless steel PEEK
	XT38253	Docklocs® Laboratory Set, up to 40° divergence compensation: 2 pcs zirconium oxide retention housing (Ø5.8 mm, height 2.5 mm) with black processing insert (height 1.9 mm), 2 pcs block-out ring, 2 pcs replacement male, green, 2 pcs replacement male, orange, 2 pcs replacement male, red	Santroprene® TPE Polyamide Housing ZrO ₂ HD-PE Purell
	XT38251	Docklocs® Laboratory Set, up to 40° divergence compensation: 2 pcs titanium retention housing (Ø5.5 mm, height 2.5 mm) with black processing insert (height 1.9 mm), 2 pcs block-out ring, 2 pcs replacement male, green, 2 pcs replacement male, orange, 2 pcs replacement male, red	Santroprene® TPE Polyamide Titanium housing HD-PE Purell

Range
















XT

	XT38205	Docklocs® Replacement Male, red Extra-light retention, 10°-20° 8 pcs	Polyamide
	XT38206	Docklocs® Replacement Male, orange Light retention, 10°-20° 8 pcs	
	XT38207	Docklocs® Replacement Male, green Strong retention, 10°-20° 8 pcs	
	XT38209	Docklocs® Block-out Ring 20 pcs	Santroprene® TPE
	XT38230	Zirconium oxide retention housing with processing insert 2 pcs	Housing ZrO ₂ HD-PE Purell
	XT38210	Titanium retention housing with processing insert 4 pcs	Titanium housing G5 HD-PE Purell
	XT38214	Docklocs® lab analogue straight (Ø4 mm) 4 pcs	Grade 5 titanium
	XT38215	Docklocs® impression coping with black processing insert 4 pcs	Titanium housing G5 HD-PE Purell

SB/RB/WB Tools	Art. No.	Name	Material
	P35601	Zeramex Rosedrill Ø2 mm	Stainless steel
	T35602	Zeradrill Pilot Ø2.3 mm	
	XT35630	Zeramex Profile Drill SB	
	T35608	Zeradrill S8 (Small 8 mm)	Stainless steel with carbon coating
	T35610	Zeradrill S10 (Small 10 mm)	
	T35612	Zeradrill S12 (Small 12 mm)	

Range

XT

	T35614	Zeradrill S14 (Small 14 mm)	Stainless steel with carbon coating	
	T35620	Zeratap Small Ø3.5 mm		
	T35622	Zeradrill Extension Small		
	T36608	Zeradrill R8 (Regular 8 mm)		
	T36610	Zeradrill R10 (Regular 10 mm)		
	T36612	Zeradrill R12 (Regular 12 mm)		
	T36614	Zeradrill R14 (Regular 14 mm)		
	T36620	Zeratap Regular Ø4.2 mm		
	T36622	Zeradrill Extension Regular		
	T37608	Zeradrill W8 (Wide 8 mm)		
	T37610	Zeradrill W10 (Wide 10 mm)		
	T37612	Zeradrill W12 (Wide 12 mm)		
	T37620	Zeratap Wide Ø5.5 mm		
	T37622	Zeradrill Extension Wide		
	C7650	Drill extension		Stainless steel

Range

XT

	T38650	Zeramex T depth gauge, 4 pcs	Stainless steel
	XT38619	Zeramex Prosthetic Key, 19 mm	
	XT38623	Zeramex Prosthetic Key, 23 mm	
	XT38628	Zeramex Prosthetic Key, 28 mm	
	XT36620	Zeramex Pick-up, 20 mm SB/RB/WB	
	XT36625	Zeramex Pick-up, 25 mm SB/RB/WB	
	XT36622	Zeramex Rescue Pick-up SB/RB/WB	
	P48932	Zeramex P Ratchet Adapter Unit Short	
	P48935	Zeramex P Surgical Ratchet without Adapter	
	XT35651	Zeramex XT Drill Stop, for Small Drill	
	XT36651	Zeramex XT Drill Stop, for Regular Drill	
	XT37651	Zeramex XT Drill Stop, for Wide Drill	
	XT48860	Zeramex XT Prosthetics Kit	
	XT48865	Zeramex XT Prosthetics Kit incl. Ratchet	

Range

XT

	XT48850	Zeramex XT Surgery Tray	
	XT48854	Zeramex XT Surgery Tray, fully equipped, incl. Ratchet	

Docklocs® is a registered trademark of MEDEALIS GmbH, DE.
Santoprene® is a registered trademark of Exxon Mobil Corporation, USA.

	0050		Caution
	Article number		Expiry date
	Batch code		Observe the instructions for use
	Do not reuse		Date of manufacture
	Non-sterile		Manufacturer
	Sterilization by steam or dry heat		Do not use if the packaging is damaged
	Keep away from direct sunlight		The products of the portfolio are FDA-approved
	Keep dry		Medical device
	Unique Device Identifier		EU authorised representative
	Sterile barrier system		Double packaging with single internal sterile barrier
	Do not resterilise		Gamma sterilization
Rx only	CAUTION: US federal law restricts the sale of this product to a physician or on his or her order.		
	Implant Small ø3.5 mm SB		Implant Regular ø4.2 mm RB
	Implant Wide ø5.5 mm WB		